

## ZN-16 Antibody

Catalog No: #62183

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	ZN-16 Antibody
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total ZN-16 protein
Immunogen Type	Peptide
Immunogen Description	The antiserum was produced against synthesized peptide derived from human ZNF292
Target Name	ZNF292
Uniprot	O60281;
Calculated MW	Predicted band size: 300 kDa
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB 1:500-2000

Note: This product is for in vitro research use only