

## RGNEF Antibody

Catalog No: #62106

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	RGNEF Antibody
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC;ICC/IF
Species Reactivity	(predicted: Human, Mouse, Rat)
Specificity	The antibody detects endogenous level of total RGNEF protein
Immunogen Type	Peptide
Immunogen Description	KLH conjugated synthetic peptide derived from human RGNEF
Uniprot	Q8N1W1;
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC 1:100-500;ICC/IF 1:100-500

Note: This product is for in vitro research use only