

## RABEP1 Antibody

Catalog No: #62192

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	RABEP1 Antibody
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC;IF
Species Reactivity	Human; Mouse; Rat
Specificity	The antibody detects endogenous level of total RABEP1 protein
Immunogen Type	Peptide
Immunogen Description	The antiserum was produced against synthesized peptide derived from human RABEP1
Target Name	RABEP1;
Uniprot	Q15276;
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC 1:100-1:300;IF 1:50-200

## Images



Formalin-fixed, paraffin-embedded human breast carcinoma tissue for RABEP1 Antibody was diluted at 1:200 dilution in immunohistochemical analysis.

Note: This product is for in vitro research use only