

## THUM1 Antibody

Catalog No: #62194

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	THUM1 Antibody
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human; Mouse
Specificity	The antibody detects endogenous level of total THUM1 protein
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from human THUM1
Target Name	THUMPD1;
Uniprot	Q9NXG2;
Calculated MW	Predicted band size: 38 kDa
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB 1:500-2000

Note: This product is for in vitro research use only