

NAT13 Antibody

Catalog No: #62201

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NAT13 Antibody
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human; Mouse
Specificity	The antibody detects endogenous level of total NAT13 protein
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from human NAT13
Target Name	NAA50; MAK3; NAT13; NAT5;
Uniprot	Q9GZZ1;
Calculated MW	Predicted band size: 19 kDa
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:500-2000

Note: This product is for in vitro research use only