

## HCP1/PCFT Antibody

Catalog No: #62203

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	HCP1/PCFT Antibody
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human; Mouse; Rat
Specificity	The antibody detects endogenous level of total HCP1/PCFT protein
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from human HCP1/PCFT
Target Name	SLC46A1; HCP1; PCFT;
Uniprot	Q96NT5;
Calculated MW	Predicted band size: 50 kDa
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB 1:500-2000

Note: This product is for in vitro research use only