

Goat Anti Mouse IgG2a Biotin Conjugated

Catalog No: #L35667



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | Goat Anti Mouse IgG2a Biotin Conjugated |
| Brief Description | Goat Polyclonal |
| Host Species | Goat |
| Clonality | Polyclonal |
| Purification | The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads. |
| Applications | WB;IHC;ICC |
| Species Reactivity | Mouse |
| Immunogen Description | Mouse IgG2a |
| Conjugates | Biotin |
| Formulation | Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

WB 1:5000 - 1:20000;IHC 1:500 - 1:1000; ICC 1:500 - 1:1000

Note: This product is for in vitro research use only