

Transmembrane protein 30A Antibody

Catalog No: #62109

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Transmembrane protein 30A Antibody
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC;IF
Species Reactivity	Mouse;Rat
Specificity	The antibody detects endogenous level of total Transmembrane protein 30A protein
Immunogen Type	Peptide
Immunogen Description	KLH conjugated synthetic peptide derived from human Transmembrane protein 30A
Uniprot	Q9NV96;
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.


Application Details

IHC 1:100-500;IF 1:100-500

Images



Formalin-fixed, paraffin-embedded mouse kidney for Transmembrane protein 30A Antibody was diluted at 1:200 dilution in immunohistochemical analysis.



Formalin-fixed, paraffin-embedded rat brain for Transmembrane protein 30A Antibody was diluted at 1:200 dilution in immunohistochemical analysis.

Note: This product is for in vitro research use only