

## HCV NS3 Antibody

Catalog No: #62112

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	HCV NS3 Antibody
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ELISA
Species Reactivity	(predicted: Hepatitis C Virus)
Specificity	The antibody detects endogenous level of total HCV NS3 protein
Immunogen Description	KLH conjugated synthetic peptide derived from HCV-NS3
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

ELISA 1:5000-10000

Note: This product is for in vitro research use only