

## EBNA 3B Antibody

Catalog No: #62114

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |   |
|-----------------------|---|
| Product Name          | EBNA 3B Antibody  |
| Brief Description     | Rabbit Polyclonal   |
| Host Species          | Rabbit  |
| Clonality             | Polyclonal  |
| Purification          | affinity purified by Protein A  |
| Applications          | WB;IHC;IF   |
| Species Reactivity    | (predicted: Epstein Barr Virus)   |
| Specificity           | The antibody detects endogenous level of total EBNA 3B protein                  |
| Immunogen Type        | Peptide   |
| Immunogen Description | KLH conjugated synthetic peptide derived from Epstein Barr Virus EBNA           |
| Calculated MW         | Predicted band size: 103 kDa  |
| Formulation           | 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.               |
| Storage               | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

## Application Details

WB 1:500-2000;IHC 1:100-500;IF 1:50-200

Note: This product is for in vitro research use only