

## EBNA 3C Antibody

Catalog No: #62115

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	EBNA 3C Antibody
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB;IHC;IF
Species Reactivity	(predicted: Human Herpesvirus 4)
Specificity	The antibody detects endogenous level of total EBNA 3C protein
Immunogen Type	Peptide
Immunogen Description	KLH conjugated synthetic peptide derived from human EBNA 3C
Calculated MW	Predicted band size: 109 kDa
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB 1:500-2000;IHC 1:100-500;IF 1:100-500

Note: This product is for in vitro research use only