

GCYA2 Antibody

Catalog No: #62208

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GCYA2 Antibody
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human; Rat
Specificity	The antibody detects endogenous level of total GCYA2 protein
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from human GCYA2
Target Name	GUCY1A2; GUC1A2; GUCSA2;
Uniprot	P33402;
Calculated MW	Predicted band size: 81 kDa
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only