

Staphylococcus aureus Antibody

Catalog No: #62117

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Staphylococcus aureus Antibody
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ELISA
Species Reactivity	(predicted: Staphylococcus aureus)
Specificity	The antibody detects endogenous level of total Staphylococcus aureus protein
Immunogen Type	protein
Immunogen Description	whole cell protein of Staphylococcus aureus Rosenbach tropina
Uniprot	P0A0L2;
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ELISA 1:5000-10000

Note: This product is for in vitro research use only