

HCMV pp65 Antibody

Catalog No: #62120

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HCMV pp65 Antibody
Brief Description	Mouse Monoclonal
Host Species	Mouse
Clonality	Monoclonal
Purification	affinity purified by Protein A
Applications	ELISA
Species Reactivity	(predicted: HCMVPP65)
Specificity	The antibody detects endogenous level of total HCMV pp65 protein
Immunogen Type	protein
Immunogen Description	Recombinant HCMV pp65 protein
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ELISA 1:5000-10000

Note: This product is for in vitro research use only