

SNAT Antibody

Catalog No: #62210

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SNAT Antibody
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC
Species Reactivity	Human; Mouse; Rat
Specificity	The antibody detects endogenous level of total SNAT protein
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from human SNAT
Target Name	AANAT; SNAT
Uniprot	Q16613;
Calculated MW	Predicted band size: 23 kDa
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:500-2000; IHC 1:50-300;

Images



Formalin-fixed, paraffin-embedded human Gastric adenocarcinoma for SNAT Antibody was diluted at 1:200 dilution in immunohistochemical analysis.

Note: This product is for in vitro research use only