

Cytochrome P450 2C11 Antibody

Catalog No: #62122

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytochrome P450 2C11 Antibody
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB;IHC;IF
Species Reactivity	Rat
Specificity	The antibody detects endogenous level of total Cytochrome P450 2C11 protein.
Immunogen Type	Peptide
Immunogen Description	KLH conjugated synthetic peptide derived from rat Cytochrome P450 2C11: 21-120/500
Conjugates	Unconjugated
Target Name	Cyp2c11
Uniprot	P08683
GeneID	29277
Calculated MW	57kDa
Concentration	1mg/ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000, IHC-P=1:100-500, IHC-F=1:100-500, IF=1:100-500, ICC/IF=1:100-500

Background

most abundant member of the p450 xenobiotic-inducible superfamily in non-stimulated liver; metabolizes xenobiotics and testosterone [RGD, Feb 2006]

Note: This product is for in vitro research use only