

PACAP-R1 Antibody

Catalog No: #62159



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PACAP-R1 Antibody
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Human; Mouse; Rat
Specificity	The antibody detects endogenous level of total PACAP-R1 protein.
Immunogen Type	Peptide
Immunogen Description	KLH-conjugated synthetic peptide of human PACAP-R1
Conjugates	Unconjugated
Target Name	ADCYAP1R1
Uniprot	P41586
GeneID	117
Calculated MW	61kDa
Concentration	1mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:500-1000

Images



All lanes : PACAP-R1 Antibody at 1/1k dilution
Lane 1 : Hela whole cell lysates
Lane 2 : SHSY5Y whole cell lysates
Lysates/proteins at 20 ug per lane.
Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 46; 51; 53; 54; 56 kDa
Observed band size: 61 kDa
Exposure time: 3 seconds

Note: This product is for in vitro research use only