

TCRB Antibody

Catalog No: #62160

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TCRB Antibody
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total TCRB protein.
Immunogen Type	Peptide
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the Central region of human TCRB. The exact sequence is proprietary.
Conjugates	Unconjugated
Target Name	TCRB
Uniprot	P04435
GeneID	6957
Calculated MW	28kDa
Concentration	1mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:500-1:1000

Images



All lanes : TCRB Antibody at 1/1k dilution
Lane 1 : Jurkat whole cell lysates
Lane 2 : MOLT4 whole cell lysates
Lysates/proteins at 20 ug per lane.
Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 13 kDa Observed band size: 28 kDa
Exposure time: 8 seconds

Note: This product is for in vitro research use only