

TRDC Antibody

Catalog No: #62161

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | TRDC Antibody |
| Brief Description | Rabbit Polyclonal |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | The antibody was purified by immunogen affinity chromatography. |
| Applications | WB;IF |
| Species Reactivity | Human |
| Specificity | The antibody detects endogenous level of total TRDC protein. |
| Immunogen Type | Peptide |
| Immunogen Description | KLH-conjugated synthetic peptide encompassing a sequence within the Central region of human TRDC. The exact sequence is proprietary. |
| Conjugates | Unconjugated |
| Target Name | TRDC |
| Uniprot | B7Z8K6 |
| GeneID | 28526 |
| Calculated MW | 27kDa |
| Concentration | 1mg/ml |
| Formulation | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide. |
| Storage | Store at 4°C in dark for 6 months |

Application Details

WB: 1:50-1:200; IF: 1:50-1:200

Note: This product is for in vitro research use only