

TRDC Antibody FITC Conjugated

Catalog No: #C62161F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TRDC Antibody FITC Conjugated
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total TRDC protein.
Immunogen Type	Peptide
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the Central region of human TRDC. The exact sequence is proprietary.
Conjugates	FITC
Target Name	TRDC
Uniprot	B7Z8K6
GeneID	28526
Calculated MW	27kDa
Concentration	1mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at 4°C in dark for 6 months

Note: This product is for in vitro research use only