

TCR alpha/beta Antibody

Catalog No: #62164



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description

Product Name	TCR alpha/beta Antibody
Brief Description	Mouse Monoclonal
Host Species	Mouse
Purification	The antibody was purified by affinity chromatography.
Applications	IF, FC
Species Reactivity	Human
Conjugates	Unconjugated
Target Name	TRAC; TRBC1
Uniprot	P01848; P01850
GeneID	28755; 28639
Formulation	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, and 0.02% sodium azide.
Storage	Store at 4°C in dark for 6 months

Note: This product is for in vitro research use only