

GPR88 Antibody

Catalog No: #62165



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	GPR88 Antibody
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total GPR88 protein.
Immunogen Type	Peptide
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the C-terminal region of human GPR88. The exact sequence is proprietary.
Conjugates	Unconjugated
Target Name	GPR88
Uniprot	Q9GZN0
GeneID	54112
Calculated MW	40kDa
Concentration	1mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:500-1:1000

Images



All lanes : GPR88 Antibody at 1/500 dilution
Lane 1 : HL60 whole cell lysates
Lysates/proteins at 20 ug per lane.
Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 40 kDa
Observed band size: 40 kDa
Exposure time: 7 seconds

Note: This product is for in vitro research use only