

TAF4 Antibody

Catalog No: #62168

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TAF4 Antibody
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB,IHC
Species Reactivity	Human, Mouse, Rat,Chicken,Pig
Specificity	The antibody detects endogenous level of total TAF4 protein.
Immunogen Type	Peptide
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human TAF4. The exact sequence is proprietary
Conjugates	Unconjugated.
Target Name	TAF4
Uniprot	O00268/Q92750
GeneID	6874; 6875
Calculated MW	110kD,91kD(Calculated)
Concentration	1mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.


Application Details

WB 1:500-1000, IHC 1:50-200

Images



All lanes : TAF4 Antibody at 1/1k dilution
 Lane 1 : A2780 whole cell lysates
 Lane 2 : SKOV3 whole cell lysates
 Lane 3 : mouse embryo lysates
 Lysates/proteins at 20 ug per lane.
 Secondary
 All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
 Predicted band size: 110;91 kDa
 Observed band size: 135;105 kDa
 Exposure time: 15 seconds



Formalin-fixed,paraffin-embedded human colon cancer and adjacent normal tissues stained for TAF4 at 1/100 dilution in immunohistochemical analysis.

Note: This product is for in vitro research use only