

B3GNT7 Antibody

Catalog No: #62221

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	B3GNT7 Antibody
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen
Applications	WB
Species Reactivity	Human, Mouse, Rat
Specificity	The antibody detects endogenous level of total B3GNT7 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from human B3GNT7
Conjugates	Unconjugated.
Target Name	B3GNT7
Uniprot	Q8NFL0
GeneID	93010
Calculated MW	44kD
Concentration	1mg/ml
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:500-2000

Note: This product is for in vitro research use only