

MCAF2 Antibody

Catalog No: #62131

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MCAF2 Antibody
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC,IF
Species Reactivity	Human,Mouse (predicted: Pig,Sheep,Cow,Dog,Horse)
Specificity	The antibody detects endogenous level of total MCAF2 protein.
Immunogen Type	Peptide
Immunogen Description	KLH conjugated synthetic peptide derived from human MCAF2: 601-682/682
Conjugates	Unconjugated.
Target Name	MCAF2
Uniprot	Q5U623
GeneID	80063
Calculated MW	76kD(Calculated)
Concentration	1mg/ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.


Application Details

IHC 1:100-500, IF 1:100-500

Images



Formalin-fixed,paraffin-embedded Mouse brain tissue stained for MCAF2 at 1/100 dilution in immunohistochemical analysis.



Formalin-fixed,paraffin-embedded Mouse brain tissue stained for MCAF2 at 1/100 dilution in immunohistochemical analysis.

Note: This product is for in vitro research use only