

Cytokeratin 84 Antibody

Catalog No: #62170

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Cytokeratin 84 Antibody
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen
Applications	WB
Species Reactivity	Human, Mouse,Rat,Monkey
Specificity	The antibody detects endogenous level of total Cytokeratin 84 protein.
Immunogen Type	Peptide
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human Cytokeratin 84. The exact sequence is proprietary
Conjugates	Unconjugated.
Target Name	Cytokeratin 84
Uniprot	Q9NSB2
GeneID	3890
Calculated MW	64kD(Calculated).
Concentration	1mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:500-1000

Images



All lanes : Cytokeratin 84 Antibody at 1/1k dilution
 Lane 1 : MCF7 whole cell lysates
 Lane 2 : H1688 whole cell lysates
 Lysates/proteins at 20 ug per lane.
 Secondary
 All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
 Predicted band size: 64 kDa
 Observed band size: 60 kDa
 Exposure time: 10 seconds

Note: This product is for in vitro research use only