

Nsp1p Antibody

Catalog No: #62132

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nsp1p Antibody
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB,IHC,IF
Species Reactivity	predicted: <i>Saccharomyces cerevisiae</i>
Specificity	The antibody detects endogenous level of total Nsp1p protein.
Immunogen Type	Peptide
Immunogen Description	KLH conjugated synthetic peptide derived from <i>Saccharomyces cerevisiae</i> Nsp1p: 451-550/823
Conjugates	Unconjugated.
Target Name	Nsp1p
Uniprot	P14907
GeneID	853409
Calculated MW	86kD
Concentration	1mg/ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:1000-2000;IHC 1:50-1:200;IF/ICC 1:100-1:500

Note: This product is for in vitro research use only