

## WBSCR14 /HRP Conjugated antibody

Catalog No: #C92920H

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	WBSCR14 /HRP Conjugated antibody
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC,IF/ICC
Species Reactivity	Rat (predicted: Human,Mouse,Pig,Cow,Chicken,Dog,Horse)
Specificity	The antibody detects endogenous level of total WBSCR14/ChREBP protein.
Immunogen Type	Protein
Immunogen Description	KLH conjugated synthetic peptide derived from human WBSCR14/ChREBP: 81-180/852
Conjugates	HRP
Target Name	WBSCR14
Uniprot	Q9NP71
GeneID	647055
Calculated MW	93kD
Concentration	1mg/ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC 1:100-500, IF/ICC 1:100-1:500

Note: This product is for in vitro research use only