

E.coli/Biotin Conjugated antibody

Catalog No: #C92921B

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	E.coli/Biotin Conjugated antibody
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ELISA
Species Reactivity	predicted: Escherichia Coli
Conjugates	Biotin
Target Name	E.coli?
Concentration	1mg/ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ELISA 1:5000-1:10000

Note: This product is for in vitro research use only