

WASH1 antibody

Catalog No: #62235



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	WASH1 antibody
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human, Mouse, Rat
Specificity	This antibody detects endogenous levels of WASH1.
Immunogen Description	Synthesized peptide derived from human WASH1 AA range: 163-213.
Target Name	WASH1
Uniprot	A8K0Z3
GeneID	100287171
Calculated MW	51kD
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:500-2000

Images



Western blot analysis of lysates from U2OS cells using WASH1 antibody at 1:1000.

Note: This product is for in vitro research use only