

Canine Glucagon Like Peptide 1 (GLP1) ELISA Kit

Catalog No: #EK20900



Package Size: #EK20900-1 48T #EK20900-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Canine Glucagon Like Peptide 1 (GLP1) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Canine

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.067ng/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:competitive inhibition

Note: This product is for in vitro research use only