

Canine Gastrin (GT) ELISA Kit

Catalog No: #EK20915



Package Size: #EK20915-1 48T #EK20915-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Canine Gastrin (GT) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Canine

Application Details

Detect Range:15.6-1000pg/mL

Sensitivity:4.28pg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:competitive inhibition

Note: This product is for in vitro research use only