

Canine follicle-stimulating hormone (FSH) ELISA Kit

Catalog No: #EK20922



Package Size: #EK20922-1 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Canine follicle-stimulating hormone (FSH) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Canine

Application Details

Detect Range: 3 mIU/mL - 210 mIU/mL

Sensitivity: 1 mIU/mL

Sample Type: Serum, Plasma, Other biological fluids

Sample Volume: 1-200 μ L

Assay Time: 1-4.5h

Detection wavelength: 450 nm

Product Description

Detection Method: competitive inhibition

Note: This product is for in vitro research use only