

# Pig Epidemic Diarrhea Virus Antibody (PEDV-Ab) ELISA Kit

Catalog No: #EK20903

Package Size: #EK20903-1 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                          |
|--------------------|----------------------------------------------------------|
| Product Name       | Pig Epidemic Diarrhea Virus Antibody (PEDV-Ab) ELISA Kit |
| Brief Description  | ELISA Kit                                                |
| Applications       | ELISA                                                    |
| Species Reactivity | Pig                                                      |

## Application Details

Detect Range:Request Information

Sensitivity:Request Information

Sample Type:Serum, Plasma, Other biological fluids

Sample Volume: 1-200 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only