

Recombinant Human Inhibin beta A chain (INHBA)

Catalog No: #AP89691



Package Size: #AP89691-1 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Inhibin beta A chain (INHBA)
Brief Description	Recombinant Protein
Host Species	Baculovirus
Immunogen Description	Human INHBA-P08476
Uniprot	P08476
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.

Note: This product is for in vitro research use only