

ZNF93 Conjugated Antibody

Catalog No: #C62262



Package Size: #C62262-AF350 100ul #C62262-AF405 100ul #C62262-AF488 100ul

#C62262-AF555 100ul #C62262-AF594 100ul #C62262-AF647 100ul

#C62262-AF680 100ul #C62262-AF750 100ul #C62262-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ZNF93 Conjugated Antibody
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of ZNF93 at Human
Immunogen Description	Synthesized peptide derived from human ZNF93 AA
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Accession No.	Swiss-Prot: P35789
Uniprot	P35789
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	68 kDa
Storage	Store at 4A C in dark for 6 months

Note: This product is for in vitro research use only