

Recombinant Thyroid peroxidase(TPO)

Catalog No: #AP89695

Package Size: #AP89695-1 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Thyroid peroxidase(TPO)
Brief Description	Recombinant Protein
Host Species	Mammalian Cell
Purification	Affinity-chromatography SDS-PAGE:β 90%
Immunogen Description	P07202 HUMAN TPO:33-846aa Source: 293 cells
Uniprot	P07202
Calculated MW	92 kDa
Tag Info	C-His
Formulation	Liquid:20 mM PBS, pH 7.4, 20% glycerol Concentration:1.2 mg/mL
Storage	Store working aliquots at 4°C for up to one week.Store at -20°C, for extended storage.

Note: This product is for in vitro research use only