

## Recombinant Vitamin D-binding protein

Catalog No: #AP89697

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Vitamin D-binding protein
Host Species	Mammalian Cell
Purification	Affinity-chromatography SDS-PAGE:β 90%
Immunogen Description	Recombinant protein with HUMAN Vitamin D-binding protein
Uniprot	P02774
Formulation	20 mM PBS, pH 7.4Liquid 2.0 mg/mL
Storage	Store working aliquots at 4°C for up to one week.Store at -20°C, for extended storage.

Note: This product is for in vitro research use only