

ILDR1 Antibody HRP Conjugated

Catalog No: #C00879H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ILDR1 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB,,IHC-P,IHC-F,ICC
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human ILDR1
Conjugates	HRP
Target Name	ILDR1
Other Names	ILDR1alpha; ILDR1beta; Immunoglobulin like domain containing receptor 1 alpha; Immunoglobulin like domain containing receptor 1 alpha; Immunoglobulin like domain containing receptor 1; Immunoglobulin like domain containing receptor 1 beta; MGC50831; ILDR1_HUMAN.
Accession No.	NCBI Gene ID286676
Uniprot	Q86SU0
GeneID	286676;
Excitation Emission	N A
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Background

ILDR1 is a putative membrane receptor. It belongs to the immunoglobulin superfamily, LISCH7 family. It is mainly expressed in prostate and to a lower extent in testis, pancreas, kidney, heart and liver.

Note: This product is for in vitro research use only