

Myelin PLP Antibody FITC Conjugated

Catalog No: #C00944F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Myelin PLP Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ICC,IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 20-70 277 derived from human Myelin proteolipid protein
Conjugates	FITC
Target Name	Myelin PLP
Other Names	PLP; PMD; HLD1; MMPL; SPG2; GPM6C; PLP DM2; Myelin proteolipid protein; Lipophilin; PLP1
Accession No.	Swiss-Prot#P60201NCBI Gene ID5354
Uniprot	P60201
GeneID	5354;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular, Cell membrane
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Background

This is the major myelin protein from the central nervous system. It plays an important role in the formation or maintenance of the multilamellar structure of myelin.

Note: This product is for in vitro research use only