

ApoA4 Antibody


Catalog No: #25293

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ApoA4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	ELISA WB
Species Reactivity	Hu C
Immunogen Type	Peptide
Immunogen Description	Raised against a 20 amino acid peptide near the carboxy terminus of chicken ApoA4.
Target Name	ApoA4
Other Names	Apolipoprotein A4, ApoA-IV
Accession No.	Swiss-Prot:P06727Gene ID:
Uniprot	P06727
Concentration	1mg/ml
Formulation	Supplied in PBS containing 0.02% sodium azide.
Storage	Can be stored at -20°C, stable for one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Images



Western blot analysis of ApoA4 in chicken small intestine tissue lysate with ApoA4 antibody at 1 ug/mL

Background

Apolipoprotein A4 (also known as ApoA-IV) is a plasma protein that is O-linked glycoprotein after proteolytic processing. ApoA4 is a major constituent of HDL and it also can be found on the surface of newly synthesized chylomicrons. The physiological function of ApoA4 is not completely identified so far, but several evidences from vitro studies indicate that this apolipoprotein involves in several steps of the reverse cholesterol transport pathway, which transports cholesterol from peripheral cells back to the liver.

Note: This product is for in vitro research use only