

Osterix Antibody HRP Conjugated

Catalog No: #C00955H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Osterix Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB,,IHC-P,IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 340-400 431 derived from human Osterix
Conjugates	HRP
Target Name	Osterix
Other Names	OSX; OI11; OI12; osterix; Transcription factor Sp7; Zinc finger protein osterix; SP7
Accession No.	Swiss-Prot#Q8TDD2NCBI Gene ID121340
Uniprot	Q8TDD2
GeneID	121340;
Excitation Emission	N A
Cell Localization	Nucleus
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Transcriptional activator essential for osteoblast differentiation. Binds to SP1 and EKLF consensus sequences and to other G C-rich sequences (By similarity).

Note: This product is for in vitro research use only