

C9a C9 Antibody HRP Conjugated

Catalog No: #C01001H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C9a C9 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB,,IHC-P,IHC-F,ICC
Species Reactivity	Hu
Immunogen Description	KLH conjugated synthetic peptide derived from human Complement component C9a
Conjugates	HRP
Target Name	C9a C9
Other Names	C9a, C9 deficiency; C9 deficiency with dermatomyositis; CO9_HUMAN; Complement component 9; Complement component 9 deficiency; Complement component C9a; Complement component C9a; C9.
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Background

C9 is a plasma protein synthesized in the liver and monocytes consisting of a single polypeptide chain. C9 is a part of the membrane attack complex (MAC), an important component of the immune system. The MAC forms upon complement system activation by invading pathogenic bacteria and consists of the four major complement proteins: C5b, C6, C7 and C8. These complement proteins bind to the outer surface of the plasma membrane of the invading cell. C9 binds to the membrane associated C5b-8 protein, which leads to the circular polymerization of 12-18 C9 molecules. These polymerized C9 molecules form a ring structure in the membrane. Molecules can then diffuse freely through this transmembrane channel, causing cell lysis and destruction of the invading bacterial cell.

Note: This product is for in vitro research use only