

PAF Antibody FITC Conjugated

Catalog No: #C01010F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PAF Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	PAF
Crossing Reactivity	PAF
Immunogen Description	KLH conjugated PAF
Conjugates	FITC
Target Name	PAF
Other Names	1-Hexadecyl-2-acetyl-sn-glycero-3-phosphocholine; 1-O-Palmityl-2-acetyl-sn-glycero-3-phosphocholine monohydrate; 2-Acetyl-1-hexadecyl-sn-glycero-3-phosphocholine; 3-sn-Phosphatidylcholine; 1-hexadecyl-2-acetyl, PAFC16.
Excitation Emission	494nm 518nm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Azelaoyl PAF is a potent PPAR γ agonist which competes for the thiazolidinedione binding site.

Note: This product is for in vitro research use only