

## AAMP Antibody FITC Conjugated

Catalog No: #C01173F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	AAMP Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt, ICC, IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 80-130 434 derived from human AAMP
Conjugates	FITC
Target Name	AAMP
Other Names	Angio-associated migratory cell protein; AAMP
Accession No.	Swiss-Prot#Q13685NCBI Gene ID14
Uniprot	Q13685
GeneID	14;
Excitation Emission	494nm 518nm
Cell Localization	Cytoplasm, Cell membrane
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

## Background

Plays a role in angiogenesis and cell migration. In smooth muscle cell migration, may act through the RhoA pathway.

Note: This product is for in vitro research use only