

Caspase-3 Antibody

Catalog No: #25442

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Caspase-3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	ELISA IHC
Species Reactivity	Hu
Specificity	At least two isoforms of Caspase-3 are known to exist.
Immunogen Type	Peptide
Immunogen Description	17 amino acid peptide near the center of human Caspase-3.
Target Name	Caspase-3
Other Names	Cysteine-aspartic acid protease 3, Caspase 3, CASP3, CPP32, CPP32B, SCA-1
Accession No.	Swiss-Prot:P42574Gene ID:836
Uniprot	P42574
GeneID	836;
Concentration	1mg/ml
Formulation	PBS containing 0.02% sodium azide.
Storage	Store at 4°C for three months and -20°C, stable for up to one year.

Application Details

Predicted MW: 30kd

Images



Immunohistochemistry of Caspase-3 in human tonsil tissue with Caspase-3 antibody at 5 ug/ml.

Note: This product is for in vitro research use only