

TNFRSF18 Antibody FITC Conjugated

Catalog No: #C01448F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | TNFRSF18 Antibody FITC Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Purified by Protein A. |
| Applications | IF |
| Species Reactivity | Hu Ms Rt |
| Immunogen Description | KLH conjugated synthetic peptide aa 185-225 228 derived from mouse TNFRSF18 |
| Conjugates | FITC |
| Target Name | TNFRSF18 |
| Other Names | AITR; Gitr; Tumor necrosis factor receptor superfamily member 18; Glucocorticoid-induced TNFR-related protein; CD357; Tnfrsf18 |
| Accession No. | Swiss-Prot#O35714NCBI Gene ID21936 |
| Uniprot | O35714 |
| GeneID | 21936; |
| Excitation Emission | 494nm 518nm |
| Cell Localization | Cytoplasm |
| Concentration | 1mg ml |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

IF=1:50-200

Background

Receptor for TNFSF18. Seems to be involved in interactions between activated T-lymphocytes and endothelial cells and in the regulation of T-cell receptor-mediated cell death. Mediated NF-kappa-B activation via the TRAF2 NIK pathway (By similarity).

Note: This product is for in vitro research use only