

## PIWI-L2 Antibody


Catalog No: #25453

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PIWI-L2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	ELISA IHC
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	14 amino acid peptide near the center of human PIWI-L2.
Target Name	PIWI-L2
Other Names	PIWI-like 2, PIWIL2, PIWIL1L, HILI, MILI
Accession No.	Swiss-Prot:Q8TC59Gene ID:55124
Uniprot	Q8TC59
GeneID	55124;
Concentration	1mg/ml
Formulation	PBS containing 0.02% sodium azide.
Storage	Store at 4°C for three months and -20°C, stable for up to one year.

## Images



Immunohistochemistry of PIWI-L2 in human brain tissue with PIWI-L2 antibody at 5 ug/mL.

Note: This product is for in vitro research use only