

ZNF521 Antibody


Catalog No: #25454

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ZNF521 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	ELISA IHC
Species Reactivity	Hu
Specificity	ZNF521 antibody is human specific.
Immunogen Type	Peptide
Immunogen Description	19 amino acid peptide near the carboxy terminus of human ZNF521.
Target Name	ZNF521
Other Names	Zinc finger protein 521, EHZF, Evi3
Accession No.	Swiss-Prot:Q96K83Gene ID:25925
Uniprot	Q96K83
GeneID	25925;
Concentration	1mg/ml
Formulation	PBS containing 0.02% sodium azide.
Storage	Store at 4°C for three months and -20°C, stable for up to one year.

Images



Immunohistochemistry of ZNF521 in human kidney tissue with ZNF521 antibody at 5 ug/ml.

Note: This product is for in vitro research use only