

IL-17 Antibody HRP Conjugated

Catalog No: #C01528H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-17 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 110-150 150 derived from rat IL-17
Conjugates	HRP
Target Name	IL-17
Other Names	IL17; CTLA8; IL-17; IL-17A; Interleukin-17A; Cytotoxic T-lymphocyte-associated antigen 8; CTLA-8; IL17A
Accession No.	Swiss-Prot#Q16552NCBI Gene ID3605
Uniprot	Q16552
GeneID	3605;
Excitation Emission	N A
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Induces stromal cells to produce proinflammatory and hematopoietic cytokines. Enhances the surface expression of ICAM1 intracellular adhesion molecule 1 in fibroblasts.

Note: This product is for in vitro research use only